



## KNEE WALKER WARNINGS & SAFETY INSTRUCTIONS



- Maximum Weight Capacity is 300 lbs. **DO NOT** exceed the maximum weight capacity (including basket contents) or serious injury could result.
- Maximum weight capacity for the Knee Walker basket is 5lbs.
- **DO NOT** use the Knee Walker on stairs, curbs, obstacles, etc. Serious risk of fall or injury may occur
- Watch out for: cords, thrown rugs, toys on the floor, spilled water or oil, or any other obstacles.
- **DO NOT** use the Knee Walker if the brake system is not working properly.
- **DO NOT** use the Knee Walker unless all parts are secure and moving parts are in good working order.
- **DO NOT** lean forward, backward, or to either side while using.
- **DO NOT** use the Knee Walker unless seat lock lever is securely locked into place.
- **DO NOT** perform any adjustments to the Knee Walker while it is being loaded into a vehicle or while it is being transported in a vehicle.
- **DO NOT** sit on the knee rest.
- **DO NOT** hang anything from the Knee Walker handles or frame at any time. This may cause the Knee Walker to tip, resulting in bodily injury or Knee Walker damage. Items should be placed in the basket.
- **DO NOT** attempt to reach objects while using the Knee Walker. Reaching for these objects will cause a change of weight distribution and may tip over, resulting in injury or damage.
- **DO NOT** walk backwards while using the Knee Walker
- ALL WHEELS MUST be in contact with the floor at all times during use. This will ensure the Knee Walker is properly balanced
- Always keep both hands on the handle grips when in use.
- When turning the Knee Walker, exercise caution and only turn while at a slow rate of movement.
- The Knee Walker is for individual use only. **DO NOT** attempt to push the Knee Walker while someone is seated on the seat. Doing so may in result in serious injury and damage to the Knee Walker.
- **The Knee Walker is not to be used as a transportation device. The Knee Walker is a walking aid only. Usage on extremely rough terrain for extended periods of time that cause undue wear on wheels and tread may cause your deposit to be withheld for payment of repairs and parts.**

**DEPOSIT NOTICE:** Deposit is refunded if Knee Walker is returned in like-condition when rented. **Extensive damage to frame, braking system, rubber grips/handles, wheels, tread and/or basket may cause the deposit to be partly or wholly withheld for repairs.** Preexisting damage will be noted on the backside of the signature slip.

I have received the **WARNINGS and SAFETY INSTRUCTIONS** sheet for **the Knee Walker** and agree to follow them. I also agree that in no event shall Spence's Pharmacy, Inc. will be liable for any direct, indirect, punitive, incidental, special consequential damages, to property or life, whatsoever arising out of or connected with the use or misuse of the rented product. I agree that any liability will be handled with the manufacturer of the product. I also agree to the pricing schedule and deposit rules.

NAME \_\_\_\_\_ DATE \_\_\_\_\_

Address \_\_\_\_\_

Phone # \_\_\_\_\_ ID \_\_\_\_\_

Deposit **\$100** Date \_\_\_\_\_

Fee Schedule:

**\$10/day \$30/week \$50/month**

**\$10** minimum charge

Date returned: \_\_\_\_\_ Initial \_\_\_\_\_

Signature \_\_\_\_\_ NOTES:

---

Notes of pre-existing damage before release:

DEPOSIT NOTICE: Deposit is refunded if Knee Walker is returned in like-condition when rented. **Extensive damage to frame, braking system, rubber grips/handles, wheels and tread and basket may cause the deposit to be partly or wholly withheld for repairs.**

initials \_\_\_\_\_

version 2018